Technical Working Party on Testing Methods and Techniques	TWM/3/20
Third Session	Original: English
Beijing, China, April 28 to May 1, 2025	Date: April 9, 2025

## THE USE OF BIOMOLECULAR TECHNOLOGY IN DUS TESTING - A CASE STUDY ON BARLEY

Document prepared by an expert from the United Kingdom

Disclaimer: this document does not represent UPOV policies or guidance

The annex to this document contains a copy of a presentation "The use of biomolecular technology in DUS testing - a case study on Barley", to be made by an expert from the United Kingdom, at the third session of the TWM.

[Annex follows]

TWM/3/20

ANNEX



















Marker selection for varietal ID and seed sample authentications Uniqueness selection with 1171 genotypes 1200 900 • Open source software "Uniqueness" **Genotypes differentiated** used to identify a subset of markers 600 based on their discrimination ability 20 SNPs can completely differentiate 300 between the 1171 genotypes in our collection 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25 Number of SNPs Each selection completed with ~ 2 million iterations of the aglorithm NIAB world-class experience, skills and resources







13

